

MODEL	ELOFLEX M-PLUS	ELOFLEX L	ELOFLEX F	ELOFLEX X	ELOFLEX H
WEIGHT <sup>1</sup>	26 kg	26 kg	26 kg	28 kg	38 kg
SPEED	6,5 km/h	6,5 km/h	6,5 km/h	6,5 km/h	8 km/h
RANGE	30 km	30 km	30 km	25 km	40 km
MOTOR	2 x Brushless motors	2 x Brushless motors	2 x Brushless motors	2 x Brushless motors	2 x Brushless motors
MOTOR EFFECT	2 x 250 Watt	2 x 200 Watt	2 x 250 Watt	2 x 250 Watt	2 x 400 Watt
BATTERY	2 x Lithium ION 10 AH / 24 V	2 x Lithium ION 10 AH / 24 V	2 x Lithium ION 10 AH / 24 V	2 x Lithium ION 7,5 AH / 24 V	2 x Lithium ION 20 AH / 24 V
BATTERY EFFECT	2 x 240 WH	2 x 240 WH	2 x 240 WH	2 x 180 WH	2 x 480 WH
BATTERY WEIGHT	1,6 kg	1,6 kg	1,6 kg	1,6 kg	3,2 kg
CHARGING TIME	3 h	3 h	3 h	3 h	5 h
TURNING RADIATION	80 cm	76 cm	78 cm	88 cm	99 cm
WHEELS	8" front – solid 12" rear – air	8" front – solid 10" rear – solid	8" front – solid 12" rear – air	8" front – solid 12" rear – air	9" front – air 13" rear – air
SEAT HIGHT	47 cm	52 cm	52 cm	51 cm	54 cm
SEAT WIDE	46 cm	46 cm	45 / 50 / 55 cm	45 / 50 / 55 cm	51 cm
SEAT DEPTH	43 cm	43 - 48 cm	43 - 48 cm	43 - 48 cm	51 cm
MAX LOAD	Max 150 kg	Max 120 kg	Max 120 kg	Max 120 kg	Max 200 kg
JOYSTICK	Left, right or backrest <sup>2</sup>	Left, right or backrest <sup>2</sup>	Left, right or backrest <sup>2</sup>	Left, right or backrest <sup>2</sup>	Left, right or backrest <sup>2</sup>
CUSHION	Removable cushion and cover	Removable cushion and cover	Removable cushion and cover	Removable cushion	Removable cushion and cover
BACKREST	Adjustable height	Adjustable trim	Adjustable angel <sup>3</sup> and trim	Adjustable angel <sup>3</sup> and trim	Adjustable trim
SIZE	90 cm (H) x 64 cm (W) x 100 cm (D)	95 cm (H) x 59 cm (W) x 99 cm (D)	95 cm (H) x 63 cm (W) x 100 cm (D)	100 cm (H) x 62 cm (W) x 102 cm (D)	103 cm (H) x 66 cm (W) x 115 cm (D)
SIZE FOLDABLE	75 cm (H) x 64 cm (W) x 40 cm (D)	70 cm (H) x 59 cm (W) x 35 cm (D)	75 cm (H) x 63 cm (W) x 40 cm (D)	76 cm (H) x 32 cm (W) x 78 cm (D)	83 cm (H) x 66 cm (W) x 45 cm (D)
GARANTIE	2 years	2 years	2 years	2 years	2 years
ARTIKELNR	7350006080029	7350006080067	7350006080685	7350006080623	7350006080679
PRICE INCL TVA	3 290 €	3 090 €	3 290 €	3 090 €	4 690 €

1) Total weight including one battery on M-Plus, L, F och H. 2) Attendant control requires holder sold as an accessory. 3) Angel 2, 12, 22 degrees.

Recommended pressure for the rear wheels = 2,5 kg