

ELOFLEX MODELS

MODEL	F	X	D2	P	H	L	S1
Weight	26 kg	28 kg	16+17 kg	28 kg	38 kg	26 kg	41 kg
Range	30 km	25 km	30 km	30 km	40 km	30 km	30 km
Max load	120 kg	120 kg	160 kg	160 kg	200 kg	120 kg	120 kg
Speed	6,5 km/h	6,5 km/h	6,5 km/h	6,5 km/h	8 km/h	6,5 km/h	9 km/h
Battery	2 x Lithium ION 10AH/24V	2 x Lithium ION 7,5AH/24V	2 x Lithium ION 10AH/24V	2 x Lithium ION 10AH/24V	2 x Lithium ION 20AH/24V	2 x Lithium ION 10AH/24V	2 x Lithium ION 10AH/24V
Turning radiation	78 cm	88 cm	89 cm	85 cm	99 cm	76 cm	200 cm
Wheels Front Rear	8" solid 12" air	8" solid 12" air	8" solid 12" air	8" solid 12" air	9" air 13" air	8" solid 10" solid	12" air 12" air
Seat hight	52 cm	51 cm	50 cm	54 cm	54 cm	52 cm	55-60 cm
Seat width	45/50/55 cm	45/50/55 cm	45/50/55 cm	45/50/55 cm	51 cm	45/50/55 cm	45/50/55 cm
Seat Depth	43-48 cm	43-48 cm	44-47 cm	43-48 cm	51 cm	43-48 cm	40 cm
Size	95x63x100 cm	100x62x102 cm	96x62x107 cm	90x62x107 cm	103x66x115 cm	95x59x99 cm	105x63x91 cm
Size folded	75x63x40 cm	76x32x78 cm	76x62x50 cm	74x62x40 cm	83x66x45 cm	70x59x35 cm	105x63x55 cm

Total weight includes one battery. Recommended pressure for the air-filled wheels = 2,5 kg. Eloflex is classified as a medical device with the Swedish Medicines Agency and meets the applicable requirements in the EU regulation (2017/745), Medical Device Regulations (MDR), which ensures the safety and performance of medical devices. It is CE approved, tested and approved according to EN12184, ISO-7176, IEC-60601 and the batteries are flight approved according to UN38.3.